

3D Printed Bands

WHY USE THESE OVER TRADITIONAL?


- no more separators
- no more guesswork
- no more long appointments

Once you've sent us your intra-oral scan, we'll be able to deliver you a highly accurate 3d printed band, custom tailored to your patient's anatomy that can be seated without the need for any separators, giving them and you, a quick and pleasant visit.

Stop worrying about bands not fitting. Stop worrying about adjustments. Stop worrying about the time.

Start getting used to appliance durability, time and cost savings, and all your very happy patients.

Just send us a scan and **we'll handle the rest**



TRADITIONAL

- 1st Appointment: Place separators
(after this appointment, patient must endure discomfort till the affected teeth have been properly spaced)
- 2nd Appointment: Take out separators, take impression, put in new separators
- 3rd Appointment: Remove separators, place fixed appliance

3D BANDS

- 1st Appointment: Take intra-oral scan of patient
(no more patient discomfort)
- 2nd Appointment: Deliver a naturally seated fixed appliance using 3d printed bands
(custom tailored fit allows easy placement that is countoured to your patient's anatomy with very little to no adjustment)

Explore more about 3D Bands:

Video for 3d bonding procedure: <https://bit.ly/3DBandsBonding>

Video

A study for 3d bands: https://bit.ly/3DBands_Article

Article